

Time Decrypt Programme for LabView Stations

Josep

May 17, 2026

Here is the guidance for Time Decrypt Programme for LabView Stations. The detail is shown in Fig.1 and Fig.2. The input is "Encrypt Injection Time" and "Key", which is upgrade in the github link everyday. The output is 3 numbers: "Injection Hour" "Injection Minute" and "Injection Second".

What you need to do is:

- Add a trigger logic in the current trigger programme you are using: At the output time of my programme, generate a 2-second square pulse.
- After confirming the programme works normally, delete the indicator of "Injection Hour" "Injection Minute" and "Injection Second", to make blind injection.




Figure 1: programme chart




Figure 2: front panel